



## Acellular dermal matrix for soft tissue reconstruction

### PRODUCT DESCRIPTION

CELLIS® Dental is a cell-free, non-pyrogenic collagen matrix (acellular dermal matrix ADM) obtained from porcine skin. CELLIS® Dental is intended for use as a surgical membrane in soft tissue repairs and serves to support, cover or replace tissue. CELLIS® Dental comes in double sterile packaging, is supplied dry and does not contain any preservatives. This surgical membrane is a resilient, biocompatible implant, which incorporates into the host tissue through cellular and microvascular infiltration and should not require another surgical procedure due to removal.

### COMPOSITION

Sterile, acellular, type I & III porcine dermis derived collagen matrix.

### INDICATIONS

CELLIS® Dental is a Class III medical device dedicated to dental surgery.

It is intended to be used in implantology, periodontology, stomatology & CMF (cranio-maxillo-facial surgery) for dental regeneration of soft tissue defects and development of soft tissues around teeth or implants.

Indications for use include:

- Treatment of recession defects (gingival recession)
- Soft tissue augmentation around teeth and implants
- Soft tissue grafting in combination with GBR/GTR
- Covering/sealing of extraction sockets

### CONTRAINDICATIONS

CELLIS® Dental must not be used in patients with:

- Known hypersensitivity or allergy to porcine materials
- Acute or chronic infection or inflammation of the buccal cavity
- General illness contra-indicated for implantology, periodontology, stomatology, cranio-maxillo-facial surgery

### REFERENCES

SIZES	REFERENCES	THICKNESSES	SIZES	REFERENCES	THICKNESSES
10 X 10 MM	NP1010E	1,4MM	20 X 50 MM	NP2050E	1,4MM
10 X 30 MM	NP1030E	1,4 MM	30 X 40 MM	NP3040E	1,4MM
15 X 20 MM	NP1520E	1,4MM	50 X 50 MM	NP5050E	1,4MM
20 X 30 MM	NP2030E	1,4MM			

## CELLIS® BENEFITS

- Quick hydration
- Easy to handle and place
- Freeze-dried
- Free from preservative



### Storage

- Store at room temperature as indicated on the label
- Keep away from heat sources and direct sunlight
- Store inside the original packaging

### Contact

#### MECELLIS BIOTECH

Address: 75 rue de Québec

17000 La Rochelle - France

Ph: +33 (0) 5 17 81 07 63

Email: [contact@meccellis.com](mailto:contact@meccellis.com)

Website: [www.meccellis.com](http://www.meccellis.com)

Brochure is intended for healthcare professionals only.  
The "Instructions for Use" should be read carefully.