



ACELLULAR DERMAL MATRIX
FOR SOFT TISSUE

RECONSTRUCTION

PRODUCT DESCRIPTION

CELLIS® Dental is a cell-free, non-pyrogenic collagen matrix (acellular dermal matrix ADM) obtained from porcine skin. CELLIS® Dental is intended for use as a surgical membrane in soft tissue repairs and serves to support, cover or replace tissue. CELLIS® Dental comes in double sterile packaging, is supplied dry and does not contain any preservatives. This surgical membrane is a resilient, biocompatible implant, which incorporates into the host tissue through cellular and microvascular infiltration and should not require another surgical procedure due to removal.



STORAGE

- Store at room temperature as indicated on the label
- Keep away from heat sources and direct sunlight
- Store inside the original packaging

LES PLUS CELLIS®




- + Quick hydration
- + Easy to handle and place

- + Freeze-dried
- + Free from preservative



**PORCINE ACELLULAR
DERMAL MATRIX FOR SOFT
TISSUE RECONSTRUCTION**



MAKING SAFE TECHNOLOGIES
AVAILABLE

REFERENCES

SIZES	REFERENCES	THICKNESSES
10 x 10 cm	NP1010E	1,4 mm
10 x 30 cm	NP1030E	1,4 mm
15 x 20 cm	NP1520E	1,4 mm
20 x 30 cm	NP2030E	1,4 mm
20 x 50 cm	NP2050E	1,4 mm
30 x 40 cm	NP3040E	1,4 mm
50 x 50 cm	NP5050E	1,4 mm

COMPOSITION

Sterile, acellular, type I & III porcine dermis derived collagen matrix.

INDICATIONS

CELLIS® Dental is a Class III medical device dedicated to dental surgery.

It is intended to be used in implantology, periodontology, stomatology & CMF (cranio-maxillo-facial surgery) for dental regeneration of soft tissue defects and development of soft tissues around teeth or implants.

Indications for use include:

- Treatment of recession defects (gingival recession)
- Soft tissue augmentation around teeth and implant
- Soft tissue grafting in combination with GBR/GTR
- Covering/sealing of extraction sockets

CONTRE-INDICATIONS

CELLIS® Dental must not be used in patients with:

- Known hypersensitivity or allergy to porcine materials
- Acute or chronic infection or inflammation of the buccal cavity
- General illness contra-indicated for implantology, periodontology, stomatology, cranio-maxillo-facial surgery

Brochure is intended for healthcare professionals only.
The « Instructions for Use » attached to the packaging should be read carefully.

Visit

www.meccellis.com

and contact your sales representative for more information.



MECCELLIS BIOTECH

75, rue de Québec
17000 La Rochelle
FRANCE

Ph. : +33 (0)5 17 81 07 63
contact@meccellis.com
www.meccellis.com



Brochure reference: TF_02_013
Version V02 - 12/11/2021