


Sub-pectoral



Pre-pectoral

ACELLULAR DERMAL MATRIX
**FOR BREAST
RECONSTRUCTION**

SURGERY



cellis[®]


Breast

TISSUE REGENERATION SUPPORTING BREAST RECONSTRUCTION

MECCCELLIS BIOTECH is an international company recognized as an expert in tissue and cell engineering.

Meccellis Biotech aims to certify high quality medical devices so that health professionals can concentrate on providing the best patient care.

MECCCELLIS BIOTECH has developed a new range of acellular dermal matrix to reinforce soft tissue for breast reconstruction. CELLIS® Breast is a surgical membrane made of porcine dermis.




Classical

CELLIS® BENEFITS






- + Quick hydration
- + Easy to handle and place
- + Designed with fenestrations to help fluid drainage
- + Adapted shapes and sizes available
- + Freeze-dried
- + Free from preservative

CELLIS® BREAST FOR SUB-PECTORAL BREAST RECONSTRUCTION

Sizes	References	Thicknesses	Shape
12 X 18 cm	CB1218F	0,9 mm	
12 X 24 cm	CB1224F	0,9 mm	







Sizes	References	Thicknesses	Shape
18 x 25 cm	CB1825FS	0,9 mm	
18 x 25 cm	CB1825FL	0,9 mm	



CELLIS® BREAST POCKET FOR PRE-PECTORAL BREAST RECONSTRUCTION

CELLIS® BREAST FOR SUB AND PRE-PECTORAL BREAST RECONSTRUCTION



Sizes	References	Thicknesses	Shape
10 x 20 cm	CB1020MF	0,9 mm	
12 x 18 cm	CB1218MF	0,9 mm	

PRODUCT DESCRIPTION

CELLIS® Breast is a cell-free, non-pyrogenic collagen matrix (acellular dermal matrix ADM) obtained from porcine skin. CELLIS® Breast is intended for use as a surgical membrane in soft tissue repairs and serves to support, cover or replace tissue. CELLIS® Breast is available in various sizes, shapes and thicknesses. CELLIS® Breast comes in double sterile packaging, is supplied dry and does not contain any preservatives. This surgical membrane is a resilient, biocompatible implant, which incorporates into the host tissue through cellular and microvascular infiltration, and should not require another surgical procedure due to removal.

COMPOSITION

Sterile, acellular, type I & III porcine dermis derived collagen matrix.

INDICATION

CELLIS® Breast is intended for implantation to reinforce soft tissue where weakness exists and for surgical repair of damaged or ruptured soft tissue membranes.

It's indicated to be used to reconstruct, to recontour and to reform the host's human soft connective tissue particularly where loss of tissue has occurred and as a supporting tissue in surgical procedures such as breast reconstruction.

CELLIS® Breast is intended for implantation for breast reconstruction after mastectomy.

CONTRAINDICATIONS

The surgical membrane is made of porcine dermis and cannot be used in patients with known hypersensitivity to porcine materials.

STORAGE

- Store at room temperature as indicated on the label.
- Keep away from heat sources and direct sunlight.
- Store inside the original packaging.

Class III medical device



date of first placing on the market: August 2018

Brochure is intended for healthcare professionals only.
The « Instructions for Use » attached to the packaging should be read carefully.

Visit

www.meccellis.com

and contact your sales representative for more information.



MECCELLIS BIOTECH

75, rue de Québec
17000 La Rochelle
FRANCE

Tél. : +33 (0)5 17 81 07 63

contact@meccellis.com
www.meccellis.com



Distributed by